

USER MANUAL

NutriMedAI

Personalized nutrition guidance
from a photograph of your meal.

Website	https://nutri-medi-ai.vercel.app
Document	User manual for people using the live application
Audience	Patients, caregivers, and wellness users

This guide explains how to sign in, set a health profile, analyze a meal, and read the result. It is written for everyday use. It does not cover installation or system administration.

SECTION 01

About this guide

NutriMedAI looks at a photo of food and explains how that meal fits the health conditions you choose. The result includes a nutrition score, estimated nutrients, short advice for your conditions, and safer alternatives.

Use this manual when you want to:

- Open the application and sign in.
- Tell the app which conditions should shape the advice.
- Upload a meal and understand each part of the result.
- Save a report or return to an earlier analysis.

Contents

- 01 About this guide
- 02 What NutriMedAI does
- 03 Before you begin
- 04 Sign in
- 05 Set your medical profile
- 06 Analyze a meal
- 07 Read your results
- 08 History, tour, and PDF reports
- 09 Practical tips
- 10 If something does not work

SECTION 02

What NutriMedAI does

The application is a web dashboard. You do not install anything. After you sign in, you describe your health context, upload a food image, and receive a written analysis tailored to that context.

Each analysis can include

- Identified food. The dish name recognized from the photograph.
- Nutrition score. A score from 0 to 100 for how the meal fits the profile you set.
- Key metrics. Estimates such as calories, protein, carbohydrates, fat, fiber, sugar, and sodium.
- Condition advice. Separate notes for a current condition and for a condition you want to monitor.
- Ingredient takeaways. A short TL;DR for the main ingredients.
- Similar but safer. Practical swaps or adjustments.

IMPORTANT

NutriMedAI provides general dietary guidance. It is not a diagnosis, a prescription, or a substitute for advice from a physician or a registered dietitian. When a result conflicts with your care plan, follow your clinician.

SECTION 03

Before you begin

You need three things:

- A current web browser on a phone or computer.
- The website address below.
- The email and password issued for your account.

Open this address	https://nutri-medi-ai.vercel.app
Email	nutriai@admin.com
Password	nutriai@admin

FIRST VISIT

If the site has been idle, the first page can take 30 to 60 seconds to respond. Wait for the page to finish loading, then sign in. Later visits in the same session are faster.

This is a shared demonstration account. Analyses saved on it are visible to anyone who signs in with the same details. Do not upload images you would not want associated with a shared login.

SECTION 04

Sign in

- 1** Open the website
Go to <https://nutri-medi-ai.vercel.app>. The home page introduces the product.
- 2** Choose Login
Select Login in the top corner. Get Started opens the same sign-in screen. Public self-registration is not available on this deployment.
- 3** Enter your details
Type the email and password, then select Log in. A successful sign-in opens the dashboard.
- 4** Leave when you are finished
Select Log out in the top bar. This is especially important on a shared computer.

If the message says the email or password is invalid, re-enter both values carefully. The password is case-sensitive. If the message says login failed, the service may still be waking up. Wait a moment and try again.

SECTION 05

Set your medical profile

The profile is what makes the advice personal. Set it before you upload a photo. You can change it for every new analysis.

Current medical condition

This is a condition you already have. Search the list or type your own. You may select more than one. Examples include diabetes, hypertension, heart disease, kidney disease, high cholesterol, celiac disease, lactose intolerance, acid reflux, thyroid disorder, and weight management. If none apply, you can leave the field unused.

Conditions to monitor

This is a concern you want the meal checked against, even if it is not a current diagnosis. Heart health and blood sugar are typical examples. Advice for this field appears under Concerned condition, separate from the current-condition advice.

HOW THE TWO FIELDS DIFFER

Current medical condition shapes advice for a condition you are already managing. Conditions to monitor shapes advice for a risk you want to watch. Use both when you need the meal reviewed from two angles.

SECTION 06

Analyze a meal

1

Start from a clear screen

If a previous result is still open, select New analysis. The profile and upload form return.

2

Confirm the profile

Check both condition fields. Add a custom condition if it is not in the list.

3

Upload the food image

Select the upload area and choose a photo. Accepted formats are PNG, JPG, JPEG, GIF, and WEBP. A preview appears after the file is chosen.

4

Add a note if it helps

The description is optional. Use it for portion size, cooking method, or a specific question, such as "one bowl, pan-fried, is the rice portion high for diabetes?"

5

Select Analyze food

The button stays inactive until an image is attached. A short progress screen appears while the meal is reviewed.

A useful photograph shows the full plate in even light, without heavy filters. Include sides and drinks if they are part of the meal. If several dishes are in the frame, say which one you want reviewed in the description.

SECTION 07

Read your results

The result opens in two areas. On the left is the identified dish, the photo, and the nutrition score. On the right is the written summary and the condition-specific advice.

ON SCREEN	WHAT IT MEANS
Identified food	The name the application assigns to the dish. If it looks wrong, run a new analysis with a clearer photo or a short description.
Nutrition score	A 0–100 reading of how the meal fits the profile on this analysis. It is a guide, not a clinical grade.
Nutrition summary	A short paragraph describing the meal in plain language.
Key metrics	Estimated calories and nutrients. Treat the numbers as estimates from a photograph, not as a laboratory assay.
Current condition	Advice written for the conditions you listed as current, including a short TL;DR and labeled points.
Concerned condition	The same style of advice for the conditions you asked the app to monitor.
Similar but safer	Swaps or portion changes that would suit the profile better.
Additional note	An extra remark when the analysis has one more point worth keeping.

Advice labels

Condition cards use short labels so you can scan them. Important marks a point to take seriously. Reasoning explains why. Action is something you can change. Benefit notes what works in your favor. Ask your doctor flags a decision that should go back to your clinician.

ESTIMATES, NOT MEASUREMENTS

A photograph cannot weigh a serving. Sauces, oil, and hidden ingredients are easy to miss. Use the metrics as a starting point, and correct them with the description field when you know the portion.

SECTION 08

History, tour, and PDF reports

Recent history

On a computer, recent analyses stay in the left column. On a phone, open them from the clock icon beside the logo. Search by dish or condition, then select an entry to open it again. Sync history refreshes the list from the account.

Take a tour

Take a tour, in the top bar, walks through the profile, the upload area, and history. You can close it at any step. It is a good first pass if someone else is using the account for the first time.

Download PDF

When a result is open, Download PDF appears in the top bar. The file is a report of that analysis, including the food image and the condition summary. Use it to keep a copy or to discuss the meal with a clinician. The download button is available only while a result is on screen.

SECTION 09

Practical tips

- Set the profile before every analysis if your question has changed. The saved history keeps the conditions used at the time.
- One plate per photo gives a cleaner identification than a crowded table.
- Name the portion when you know it. "Two rotis and a small bowl of dal" is more useful than a photo alone.
- Review Current condition and Concerned condition separately. They answer different questions.
- Download the PDF before you start a new analysis if you want a file of the current result.
- Sign out on a shared device when you finish.

SECTION 10

If something does not work

ON SCREEN	WHAT IT MEANS
The site will not open	Wait up to a minute and refresh. The service pauses after inactivity and needs a moment to start.
Login failed	The application could not reach the sign-in service. Wait, refresh, and try once more.
Invalid email or password	Check the email and password in Section 03. The password is case-sensitive.
Analyze food stays disabled	An image has not been attached yet. Choose a file and wait for the preview.
The dish name is wrong	Upload a closer photo or name the dish in the description, then analyze again.
The result looks incomplete	Run the analysis again. A weak connection can interrupt the response.
No Download PDF button	Open a completed result first. The button is hidden on the empty upload screen.

MEDICAL NOTICE

Do not change medication, insulin, or a prescribed diet because of a score in this application. NutriMedAI supports everyday food decisions. Clinical decisions belong with your care team.